

Quick reference guide


Selecting images

Select images by holding down the Ctrl key and clicking the left mouse button with the mouse cursor over the image to be selected.













Many image management features can be accessed by clicking the right mouse button with the mouse cursor placed over the image.

Selecting between record stages

Select between record stages by clicking on the down arrow on the right side of the record stage box.



Main screen - Horizontal button bar

- | | |
|--|---|
|  Clear patient details |  Edit record details |
|  Add new patient file |  Add notes to patient's record |
|  Find/open patient file |  Cut |
|  Edit patient details |  Copy |
|  Add patient's keywords |  Paste |
|  Add new record |  Save |

Main screen - Vertical button bar















- | | |
|--|---|
|  Add images |  Add image notes |
|  Delete all images |  Play macro 1 |
|  Edit images |  Play macro 2 |
|  Image information |  Play macro 3 |
|  View full screen |  Restore to original image |
|  Create/view montages |  Revert to edited image |
|  Add image keywords |  Export selected images |

Image editing screen

- | | |
|--|---|
|  Zoom in / zoom 1:1 |  Image palette / Red-eye |
|  Rotate 90° CCW / CW |  Previous / Next image |
|  Free rotate / Crop |  Undo last/all image edits |
|  Flip horizontal / vertical |  Record macros |